
Human IL-6 Mouse mAb, PE conjugated

Cat. No:	V0406-PE
Name:	Human IL-6 Mouse mAb, PE conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6G9
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA
	Excitation spectrum: 496/564nm
	Emission spectrum: 578nm
Antigen:	Recombinant Human IL-6 Protein