

bs-2247R

[Primary Antibody]

AHA3 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>Target: AHA3</p> <p>Immunogen: KLH conjugated synthetic peptide derived from Arabidopsis thaliana (thale cress) AHA3: 851-949/949.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p>	<p>Isotype: IgG</p> <p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Arabidopsis Thaliana)</p> <p>Predicted MW: 104 kDa</p>
--	--

— SELECTED CITATIONS —

- **[IF=11.6]** Tian Ye. et al. The zinc finger protein DHHC09 S-acylates the kinase STRK1 to regulate H2O2 homeostasis and promote salt tolerance in rice. PLANT CELL. 2024 Jan;; WB ;Rice. 38180963