

bs-2247R

[Primary Antibody]

AHA3 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500)
Clonality: Polyclonal		IHC-F (1:100-500)
Target: AHA3		IF (1:100-500)
Immunogen: KLH conjugated synthetic peptide derived from Arabidopsis thaliana (thale cress) AHA3: 851-949/949.		ELISA (1:5000-10000)
Purification: affinity purified by Protein A		Reactivity: (predicted: Arabidopsis Thaliana)
Concentration: 1mg/ml		Predicted MW.: 104 kDa
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=11.6]** Tian Ye. et al. The zinc finger protein DHHC09 S-acylates the kinase STRK1 to regulate H2O2 homeostasis and promote salt tolerance in rice. PLANT CELL. 2024 Jan;; WB ;Rice. 38180963