

Human PGI Mouse mAb

Cat. No: V4101
Name: Human PGI Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 4A8
Host/Source: Mouse
Subtype: IgG1
Reactivity: Human
Purity: >90% as determined by SDS-PAGE
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA, ICA
Antigen: Recombinant Human PGA4 Protein ([bs-41228P](#))
Pairing Information:

| ELISA | | Detection | |
|---------|-------|-----------|-------|
| | | V4101 | V4102 |
| Capture | V4101 | - | +++ |
| | V4102 | - | - |
| ICA | | Detection | |
| | | V4101 | V4102 |
| Capture | V4101 | - | - |
| | V4102 | +++ | - |