

## Human PGI Mouse mAb

Cat. No: V4102  
Name: Human PGI Mouse mAb  
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  
Clone type: Monoclonal  
Clone No: 6B7  
Host/Source: Mouse  
Subtype: IgG1  
Reactivity: Human  
Format: Liquid  
Buffer: 0.01M PBS (pH7.4).  
Applications: ELISA, ICA  
Antigen: Recombinant Human PGA4 Protein ([bs-41228P](#))  
Pairing Information:

| ELISA   |       | Detection |       |
|---------|-------|-----------|-------|
|         |       | V4101     | V4102 |
| Capture | V4101 | -         | +++   |
|         | V4102 | -         | -     |
| ICA     |       | Detection |       |
|         |       | V4101     | V4102 |
| Capture | V4101 | -         | -     |
|         | V4102 | +++       | -     |