

Human PGII Mouse mAb

Cat. No: V4204
Name: Human PGII Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 5C2
Host/Source: Mouse
Subtype: IgG1
Reactivity: Human
Purity: >90% as determined by SDS-PAGE
Format: Liquid
Buffer: PBS (pH7.4).
Applications: ELISA
Antigen: Recombinant Human PGC Protein ([bs-43143P](#))
Pairing Information:

ELISA		Detection				
		V4201	V4202	V4203	V4204	V4205
Capture	V4201	-	+++	-	-	-
	V4202	-	-	-	-	-
	V4203	-	+++	-	-	-
	V4204	-	-	-	-	+++
	V4205	-	-	-	-	-