

## Human PGII Mouse mAb

Cat. No: V4205  
Name: Human PGII Mouse mAb  
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  
Clone type: Monoclonal  
Clone No: 7D5  
Host/Source: Mouse  
Subtype: IgG1  
Reactivity: Human  
Purity: >90% as determined by SDS-PAGE  
Format: Liquid  
Buffer: PBS (pH7.4).  
Applications: ELISA  
Antigen: Recombinant Human PGC Protein ([bs-43143P](#))  
Pairing Information:

ELISA		Detection				
		V4201	V4202	V4203	V4204	V4205
Capture	V4201	-	+++	-	-	-
	V4202	-	-	-	-	-
	V4203	-	+++	-	-	-
	V4204	-	-	-	-	+++
	V4205	-	-	-	-	-