

Human sST2 Mouse mAb

Cat. No:	V3705
Name:	Human sST2 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	5H4
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	Recombinant Human sST2 Protein (bs-43005P)
Pairing Information:	

ELISA		Detection	
		V3704	V3705
Capture	V3704	-	+++
	V3705	-	-
ICA		Detection	
		V3704	V3705
Capture	V3704	-	+++
	V3705	-	-