

## Human Renin Mouse mAb

Cat. No: V4301  
 Name: Human Renin Mouse mAb  
 Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  
 Clone type: Monoclonal  
 Clone No: 3F5  
 Host/Source: Mouse  
 Subtype: IgG1  
 Reactivity: Human  
 Purity: >90% as determined by SDS-PAGE  
 Format: Liquid  
 Buffer: PBS (pH7.4).  
 Applications: ELISA  
 Antigen: Recombinant Human Renin Protein ([bs-43146P](#))  
 Pairing Information:

ELISA		Detection	
		V4301	V4302
Capture	V4301	-	+++
	V4302	-	-