

## Recombinant Bovine Pregnancy-associated glycoprotein 1 (HEK293)

Catalog Number: bs-43189P

Concentration: >1mg/ml

Species: Bovine

AA Seq: 16-380/380

Predicted MW: 43 kDa

Tags: C-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE


Purification: AC

Form: Lyophilized or Liquid


Storage: PBS (pH7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

### VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.