

bs-15242R

[Primary Antibody]

C6orf223 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: ELISA (1:5000-10000)
Clonality: Polyclonal		Reactivity: (predicted: Human)
GeneID: 221416	SWISS: Q8N319	
Target: C6orf223		
Immunogen: KLH conjugated synthetic peptide derived from human c6orf223: 1-100/242.		Predicted MW.: 26 kDa
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		