

bs-4361R

[Primary Antibody]

ZNF695 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 57116</p> <p>Target: ZNF695</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human ZNF695: 101-200/515.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: ZNF695 (Zinc Finger Protein 695) is a Protein Coding gene. GO annotations related to this gene include nucleic acid binding. An important paralog of this gene is ZNF98.</p>	<p>Isotype: IgG</p> <p>SWISS: Q8IW36</p>	<p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human)</p> <p>Predicted MW.: 54 kDa</p> <p>Subcellular Location: Nucleus</p>
--	--	---