

---

## ZNF695 Rabbit pAb

Catalog Number: bs-4361R

Target Protein: ZNF695

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),  
ELISA (1:5000-10000)

Reactivity: (predicted:Human)

Predicted MW: 54 kDa

Entrez Gene: 57116

Swiss Prot: Q8IW36

Source: KLH conjugated synthetic peptide derived from human ZNF695: 101-200/515.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: ZNF695 (Zinc Finger Protein 695) is a Protein Coding gene. GO annotations related to this gene include nucleic acid binding. An important paralog of this gene is ZNF98.