
ZNF85 Rabbit pAb

Catalog Number: bs-4366R

Target Protein: ZNF85

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500)

Reactivity: (predicted:Human)

Predicted MW: 69 kDa

Entrez Gene: 7639

Swiss Prot: Q03923

Source: KLH conjugated synthetic peptide derived from human ZNF85: 1-100/595.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: ZNF85 (Zinc Finger Protein 85) is a Protein Coding gene. GO annotations related to this gene include nucleic acid binding and transcription corepressor activity. An important paralog of this gene is ZNF724.