

bs-15319R

[Primary Antibody]

C9orf171 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal GeneID: 389799 Target: C9orf171 Immunogen: KLH conjugated synthetic peptide derived from human C9orf171 : 201-300/320. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: CFAP77 (Cilia And Flagella Associated Protein 77) is a Protein Coding gene.	Isotype: IgG SWISS: Q6ZQR2	Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) Reactivity: (predicted: Human, Mouse, Rat, Sheep, Cow, Horse) Predicted MW.: 36 kDa Subcellular Location: Cytoplasm ,Nucleus
---	--	---