

---

## ZSCAN1 Rabbit pAb

Catalog Number: bs-4384R

Target Protein: ZSCAN1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human)

Predicted MW: 46 kDa

Entrez Gene: 284312

Swiss Prot: Q8NBB4

Source: KLH conjugated synthetic peptide derived from human ZSCAN1: 211-310/408.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Background:** ZSCAN1 (Zinc Finger And SCAN Domain Containing 1) is a Protein Coding gene. Diseases associated with ZSCAN1 include Arrhythmogenic Right Ventricular Dysplasia 1. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and RNA polymerase II transcription factor activity, sequence-specific DNA binding. An important paralog of this gene is ZNF274.