

bs-4384R**[Primary Antibody]****ZSCAN1 Rabbit pAb**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human) Predicted MW.: 46 kDa Subcellular Location: Nucleus
Clonality: Polyclonal		
GeneID: 284312	SWISS: Q8NBB4	
Target: ZSCAN1		
Immunogen: KLH conjugated synthetic peptide derived from human ZSCAN1: 211-310/408.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: ZSCAN1 (Zinc Finger And SCAN Domain Containing 1) is a Protein Coding gene. Diseases associated with ZSCAN1 include Arrhythmogenic Right Ventricular Dysplasia 1. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and RNA polymerase II transcription factor activity, sequence-specific DNA binding. An important paralog of this gene is ZNF274.		