
Recombinant human LPPR2 protein, N-Trx&His

Catalog Number: bs-42291P

Concentration: >0.5mg/ml

Species: Human

AA Seq: 287-343/343

Predicted MW: 24.2

Tags: N-Trx&His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: affinity purified by Protein A

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0)

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.