
Recombinant mouse Slco1a1 protein, N-Trx&His

Catalog Number: bs-42267P

Concentration: >0.5mg/ml

Species: Mouse

AA Seq: 108-158/670

Predicted MW: 23.8

Tags: N-Trx-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: mediates the transport of sulfated, amidated bile acid, sulfolithocholytaurine, into rat hepatocytes

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.