

Cardiac Troponin I (cTnI) Mouse mAb

Cat. No:	V3201
Name:	Cardiac Troponin I (cTnI) Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	17C8L
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted: Human)
Format:	Liquid
Concentration:	>1mg/ml
Buffer:	0.01M PBS (pH7.4) with 0.02% Proclin300.
Applications:	ELISA
Antigen:	Recombinant human Cardiac Troponin I (cTnI) (bs-41212P)
Pairing Information:	

ELISA		Detection				
		V3201	V3401	V3501	V3201 + V3401	V3401 + V3501
Capture	V3201	-	-	+++	-	-
	V3401	+++	-	-	-	-
	V3501	-	-	-	-	-
	V3201 + V3401	-	-	+++	-	-
	V3401 + V3501	+++	-	-	-	-