

Human HE4 Mouse mAb

Cat. No:	V7503
Name:	Human HE4 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	5D12
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Concentration:	>1mg/ml
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	Recombinant Human WAP four-disulfide core domain protein 2 Protein (bs-43151P)
Pairing Information:	

ELISA		Detection				
		V7501	V7502	V7503	V7504	V7505
Capture	V7501	-	+++	-	-	-
	V7502	-	-	-	-	-
	V7503	-	+++	-	-	-
	V7504	-	-	-	-	-
	V7505	-	-	-	-	-
ICA		Detection				
		V7501	V7502	V7503	V7504	V7505
Capture	V7501	-	-	-	-	-
	V7502	-	-	-	-	-
	V7503	-	-	-	-	-
	V7504	-	-	-	-	-
	V7505	-	-	+++	+++	-