

bs-40296G-FITC

[Secondary Antibodies]

AffiniPure Goat Anti-Mouse IgG H&L(min X Hu, Bov, Hrs, Rb, Sw Sr Prot), FITC conjugated



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— **DATASHEET** —

Host: Goat

Isotype: IgG

Clonality: Polyclonal

Target: AffiniPure Goat Anti-Mouse IgG H&L(min X Hu, Bov, Hrs, Rb, Sw Sr Prot)

Purification: affinity purified by Protein G

Storage: Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.

Applications: IF (1:100-1000)

Flow-Cyt (1:100-1000)

ICC/IF (1:100-1000)

Excitation Spectrum: 495nm

Emission spectrum: 519nm

Reactivity: Mouse