

---

## AffiniPure Goat Anti-Mouse IgG H&L(min X Hu, Bov, Hrs, Rb, Sw Sr Prot), FITC conjugated

Catalog Number: bs-40296G-FITC

Target Protein: AffiniPure Goat Anti-Mouse IgG H&L(min X Hu, Bov, Hrs, Rb, Sw Sr Prot)

Host: Goat

Clonality: Polyclonal

Isotype: IgG

Applications: Excitation spectrum: 495nm

Emission spectrum: 519nm

Not yet tested in other applications.

Optimal working dilutions must be determined by the end user.

Reactivity: Mouse

Purification: affinity purified by Protein G

Storage: Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.