

bs-40296G-FITC

[Secondary Antibodies]



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

AffiniPure Goat Anti-Mouse IgG H&L(min X Hu, Bov, Hrs, Rb, Sw Sr Prot), FITC conjugated

DATASHEET

Host: Goat

Isotype: IgG

Reactivity: Mouse

Clonality: Polyclonal

Target: AffiniPure Goat Anti-Mouse IgG H&L(min X Hu, Bov, Hrs, Rb, Sw Sr Prot)

Purification: affinity purified by Protein G

Storage: Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.