
Cystatin-C (Capture Antibody) Mouse mAb

Cat. No:	V0501
Name:	Cystatin-C (Capture Antibody) Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	1G3K
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Format:	Liquid
Concentration:	≥1mg/ml
Buffer:	0.01M PBS (pH7.4) with 0.02% Proclin300.
Applications:	ELISA
Antigen:	Recombinant human Cystatin-C (bs-41175P)