
NGAL/Lipocalin 2 Mouse mAb, Biotin conjugated

Cat. No:	V0603-Bio
Name:	NGAL/Lipocalin 2 Mouse mAb, Biotin conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3A1
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted: Human)
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Applications:	ELISA, ICA
Antigen:	Recombinant human NGAL protein