

---

## Recombinant mouse FCGR1A/CD64 protein, (HEK293)

Catalog Number: bs-43511P

Concentration: >0.5mg/ml

Species: Mouse

AA Seq: 25-297/404

Predicted MW: 34

Tags: C-His-Avi

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: PBS (pH=7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

**Background:** This gene encodes a protein that plays an important role in the immune response. This protein is a high-affinity Fc-gamma receptor. The gene is one of three related gene family members located on chromosome 1. [provided by RefSeq, Jul 2008]

### VALIDATION IMAGES

---



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.