

---

## Recombinant human GIMAP6 protein, N-His

Catalog Number: bs-42315P

Concentration: >1mg/ml

Species: Human

AA Seq: 1-305/305

Predicted MW: 36.4 kDa

Tags: N-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: affinity purified by Protein A

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH=8.0) with 8M Urea.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Highly expressed in spleen, lymph nodes, lung and placenta. Expressed at moderate level in thymus, kidney, heart and digestive tract. Weakly expressed in other lymphoid tissues.

### VALIDATION IMAGES

---



The purity of the protein is greater than 85% as determined by reducing SDS-PAGE.