
Recombinant Mumps virus Non-structural protein V, N-Trx&His

Catalog Number: bs-42330P

Concentration: >1mg/ml

Species: Mumps virus

AA Seq: 50-224/224

Predicted MW: 24

Tags: N-Trx&His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH=8.0).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

VALIDATION IMAGES



The purity of the protein is greater than 85% as determined by reducing SDS-PAGE.