

---

## Recombinant Mumps virus Fusion glycoprotein F0 protein, N-His

Catalog Number: bs-42329P

Concentration: >1mg/ml

Species: Mumps virus

AA Seq: 20-482/538

Predicted MW: 52.8

Tags: N-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH=8.0) with 8M Urea.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

### VALIDATION IMAGES

---



The purity of the protein is greater than 85% as determined by reducing SDS-PAGE.