

---

## Human VEGF165 Mouse mAb, HRP conjugated

Cat. No:	V2604-HRP
Name:	Human VEGF165 Mouse mAb, HRP conjugated
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	5D7
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted:&nbsp;Human)
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Antigen:	Recombinant human VEGF165 Protein