

## Human NT-proBNP Mouse mAb

Cat. No:	V3309
Name:	Human NT-proBNP Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	4C6
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Lyophilized
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA
Antigen:	KLH conjugated synthetic peptide derived from human NPPB ( <a href="#">bs-43095P</a> )
Pairing Information:	

ELISA		Detection					
		V3305	V3306	V3307	V3308	V3309	V3310
Capture	V3305	-	+++	-	-	-	-
	V3306	-	-	-	-	-	-
	V3307	-	-	-	-	-	-
	V3308	-	-	+++	-	-	-
	V3309	+	-	-	-	-	-
	V3310	+	-	-	-	-	-