

Human NT-proBNP Mouse mAb

Cat. No: V3310
Name: Human NT-proBNP Mouse mAb
Storage: Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 10G6
Host/Source: Mouse
Subtype: IgG1
Reactivity: Human
Purity: >90% as determined by SDS-PAGE
Format: Liquid
Buffer: 0.01M PBS(pH7.4).
Applications: ELISA
Pairing Information:

ELISA		Detection					
		V3305	V3306	V3307	V3308	V3309	V3310
Capture	V3305	-	+++	-	-	-	-
	V3306	-	-	-	-	-	-
	V3307	-	-	-	-	-	-
	V3308	-	-	+++	-	-	-
	V3309	+	-	-	-	-	-
	V3310	+	-	-	-	-	-