

---

## Human NT-proBNP Mouse mAb

Cat. No:	V3314
Name:	Human NT-proBNP Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3F12
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA
Antigen:	KLH conjugated synthetic peptide derived from human NPPB ( <a href="#">bs-43095P</a> )