
Human BNP Mouse mAb

Cat. No:	V3324
Name:	Human BNP Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	10A9
Host/Source:	Mouse
Subtype:	IgG2b
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4).
Antigen:	KLH conjugated synthetic peptide derived from human BNP Protein (bs-43095P)