

## Human BNP Mouse mAb

Cat. No:	V3325
Name:	Human BNP Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	5H8
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA
Antigen:	KLH conjugated synthetic peptide derived from human BNP Protein ( <a href="#">bs-43095P</a> )
Pairing Information:	

ELISA		Detection			
		V3321	V3322	V3325	V3326
Capture	V3321	-	+++	-	-
	V3322	-	-	-	-
	V3325	+	-	-	-
	V3326	+	-	-	-