

## Human IGFBP3 Mouse mAb

Cat. No:	V7303
Name:	Human IGFBP3 Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3B6
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human IGFBP3 Protein ( <a href="#">bs-43124P</a> )
Pairing Information:	

ELISA		Detection			
		V7301	V7302	V7303	V7304
Capture	V7301	-	+++	-	-
	V7302	-	-	-	-
	V7303	-	-	-	+++
	V7304	-	-	-	-