

## Human IGFBP3 Mouse mAb

Cat. No: V7305  
Name: Human IGFBP3 Mouse mAb  
Storage: Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.  
Clone type: Monoclonal  
Clone No: 1E9  
Host/Source: Mouse  
Subtype: IgG2b  
Reactivity: Human  
Purity: >90% as determined by SDS-PAGE  
Format: Liquid  
Buffer: 0.01M PBS(pH7.4).  
Applications: ICA  
Antigen: Recombinant Human IGFBP3 Protein ([bs-43124P](#))  
Pairing Information:

| ICA     |       | Detection |       |       |       |
|---------|-------|-----------|-------|-------|-------|
|         |       | V7305     | V7306 | V7307 | V7308 |
| Capture | V7305 | -         | -     | -     | -     |
|         | V7306 | -         | -     | -     | +++   |
|         | V7307 | +++       | -     | -     | -     |
|         | V7308 | -         | -     | -     | -     |