

Human IGFBP3 Mouse mAb

Cat. No: V7307
Name: Human IGFBP3 Mouse mAb
Storage: Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 5F6
Host/Source: Mouse
Subtype: IgG1
Reactivity: Human
Purity: >90% as determined by SDS-PAGE
Format: Liquid
Buffer: 0.01M PBS(pH7.4).
Applications: ICA
Antigen: Recombinant Human IGFBP3 Protein ([bs-43124P](#))
Pairing Information:

ICA		Detection			
		V7305	V7306	V7307	V7308
Capture	V7305	-	-	-	-
	V7306	-	-	-	+++
	V7307	+++	-	-	-
	V7308	-	-	-	-