

Human IGFBP3 Mouse mAb

Cat. No: V7308
 Name: Human IGFBP3 Mouse mAb
 Storage: Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
 Clone type: Monoclonal
 Clone No: 5H6
 Host/Source: Mouse
 Subtype: IgG2a
 Reactivity: Human
 Purity: >90% as determined by SDS-PAGE
 Format: Liquid
 Buffer: 0.01M PBS(pH7.4).
 Applications: ICA
 Antigen: Recombinant Human IGFBP3 Protein ([bs-43124P](#))
 Pairing Information:

ICA		Detection			
		V7305	V7306	V7307	V7308
Capture	V7305	-	-	-	-
	V7306	-	-	-	+++
	V7307	+++	-	-	-
	V7308	-	-	-	-