

---

## Alpha-synuclein Mouse mAb

Cat. No:	V5701
Name:	Alpha-synuclein Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	2A12
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4).
Applications:	ELISA
Antigen:	Recombinant human Alpha-synuclein protein ( <a href="#">bs-42005P</a> )