

Human NT-proBNP Mouse mAb

Cat. No:	V3308
Name:	Human NT-proBNP Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	5B2
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human ProBNP Protein (bs-43095P)
Pairing Information:	

ELISA		Detection					
		V3305	V3306	V3307	V3308	V3309	V3310
Capture	V3305	-	+++	-	-	-	-
	V3306	-	-	-	-	-	-
	V3307	-	-	-	-	-	-
	V3308	-	-	+++	-	-	-
	V3309	+	-	-	-	-	-
	V3310	+	-	-	-	-	-