

---

## **Cardiac Troponin I-C (cTnI-C) Mouse mAb, Biotin conjugated**

Cat. No:	V3501-Bio
Name:	Cardiac Troponin I-C (cTnI-C) Mouse mAb, Biotin conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	5D12K
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Concentration:	1mg/ml
Buffer:	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Antigen:	Recombinant human Cardiac Troponin I-C