
WWL113

产品编号: D50828

CAS: 947669-86-5

分子式: C₂₉H₂₆N₂O₄

纯度: ≥98%

InChi: InChI=1S/C₂₉H₂₆N₂O₄/c1-3-34-28(32)25-9-7-22(8-10-25)23-11-13-27(14-12-23)35-29(33)31(2)
20-21-5-4-6-26(19-21)24-15-17-30-18-16-24/h4-19H,3,20H₂,1-2H₃

InChi Key: AKIIPHDGVCFVCC-UHFFFAOYSA-N

Smiles: CCOC(=O)C1=CC=C(C=C1)C1C=CC(=CC=1)OC(=O)N(C)CC1C=CC=C(C=1)C1=CC=NC=C1

外观: 固体粉末

作用通路: Others

溶解性: DMSO up to 50 mM

保存条件: Store in dry, dark place for one year.

产品介绍: WWL113 is a novel potent and selective inhibitor of Ces3 and Ces1f with IC₅₀ ~0.1 μM in cellular assays.