
SMER3

产品编号: D50911

CAS: 67200-34-4

分子式: C₁₁H₄N₄O₂

纯度: ≥98%

InChi: InChI=1S/C₁₁H₄N₄O₂/c16-9-6-4-2-1-3-5(6)7-8(9)13-11-10(12-7)14-17-15-11/h1-4H

InChi Key: SFSSAKVWCKFRHE-UHFFFAOYSA-N

Smiles: O=C1C2=NC3=NON=C3N=C2C2C=CC=CC=21

外观: 固体粉末

作用通路: E1/E2/E3 Enzyme

溶解性: DMSO up to 50 mM

保存条件: Store in dry, dark place for one year.

产品介绍: SMER3 is a small molecule enhancer of rapamycin (SMER). It is also a specific inhibitor of Ubiquitin E3 ligase.