
Go6893

产品编号: D50931

CAS: 133053-19-7

分子式: C₂₆H₂₆N₄O₃

纯度: ≥98%

InChi: InChI=1S/C₂₆H₂₆N₄O₃/c1-29(2)11-6-12-30-15-20(18-13-16(33-3)9-10-22(18)30)24-23(25(31)28-26(24)32)19-14-27-21-8-5-4-7-17(19)21/h4-5,7-10,13-15,27H,6,11-12H₂,1-3H₃, (H,28,31,32)

InChi Key: LLJJDLHGZUOMQP-UHFFFAOYSA-N

Smiles: CN(C)CCCN1C=C(C2=C(C3=CNC4=CC=CC=C34)C(=O)NC2=O)C2=CC(=CC=C12)OC

外观: 固体粉末

作用通路: PKC

溶解性: DMSO up to 100 mM

保存条件: Store in dry, dark place for one year.

产品介绍: Go6893 is a potent pan-PKC inhibitor against for PKC α , PKC β , PKC γ , PKC δ and PKC ζ with IC₅₀ of 7 nM, 7 nM, 6 nM, 10 nM and 60 nM, respectively. It is inactive to PKC μ .