
CPI-169 racemate

产品编号: D50933

CAS: 1450655-76-1

分子式: C₂₇H₃₆N₄O₅S

纯度: ≥98%

InChi: InChI=1S/C₂₇H₃₆N₄O₅S/c1-6-37(34,35)30-13-11-20(12-14-30)18(3)31-19(4)25(21-9-7-8-10-23(21)31)27(33)28-16-22-24(36-5)15-17(2)29-26(22)32/h7-10,15,18,20H,6,11-14,16H₂,1-5H₃, (H,28,33)(H,29,32)

InChi Key: LHGUZCKPFXXVPV-UHFFFAOYSA-N

Smiles: CC(C1CCN(CC1)S(=O)(=O)CC)N1C2=CC=CC=C2C(C(=O)NCC2=C(C=C(C)NC2=O)OC)=C1C

外观: 固体粉末

作用通路: Epigenetic Reader Domain

溶解性: DMSO up to 100 mM, EtOH up to 100 mM

保存条件: Store in dry, dark place for one year.

产品介绍: CPI-169 is a highly potent and selective small molecule inhibitor of histone methyltransferase EZH2.