
Tenovin-3

产品编号: D51036

CAS: 1011301-27-1

分子式: C₁₈H₂₁N₃OS

纯度: ≥98%

InChi: InChI=1S/C₁₈H₂₁N₃OS/c1-18(2,3)13-6-4-12(5-7-13)16(22)21-17(23)20-15-10-8-14(19)9-11-15/h4-11H,19H2,1-3H3,(H2,20,21,22,23)

InChi Key: NLAXTZPUTNGRDU-UHFFFAOYSA-N

Smiles: CC(C)(C)C1C=CC(=CC=1)C(=O)NC(=S)NC1C=CC(N)=CC=1

外观: 固体粉末

作用通路: MDM-2/p53

溶解性: DMSO: ≥31 mg/mL

保存条件: Store in dry, dark place for one year.

产品介绍: Tenovin-3 is a sirtuin inhibitor.