

bs-4390R

[Primary Antibody]

ZMAT5 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 55954</p> <p>Target: ZMAT5</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human ZMAT5: 51-150/170.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: ZMAT5 (Zinc Finger Matrin-Type 5) is a Protein Coding gene. Among its related pathways are Gene Expression and mRNA Splicing - Major Pathway.</p>	<p>Isotype: IgG</p> <p>SWISS: Q9UDW3</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Dog)</p> <p>Predicted MW.: 20 kDa</p> <p>Subcellular Location: Nucleus</p>
--	--	--