
ZMAT5 Rabbit pAb

Catalog Number: bs-4390R

Target Protein: ZMAT5

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Pig, Dog)

Predicted MW: 20 kDa

Entrez Gene: 55954

Swiss Prot: Q9UDW3

Source: KLH conjugated synthetic peptide derived from human ZMAT5: 51-150/170.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: ZMAT5 (Zinc Finger Matrin-Type 5) is a Protein Coding gene. Among its related pathways are Gene Expression and mRNA Splicing - Major Pathway.